



# Under-the-Counter Cabinet

## Installation and Maintenance Instructions

### Cabinet Model Numbers:

Model No. 1424	Partially Lined Acid Cabinet
Model No. 1430	Partially Lined Acid Cabinet
Model No. 1436	Partially Lined Acid Cabinet
Model No. 1448	Partially Lined Acid Cabinet
Model No. 1624	Fully Lined Acid Cabinet
Model No. 1630	Fully Lined Acid Cabinet
Model No. 1636	Fully Lined Acid Cabinet
Model No. 1648	Fully Lined Acid Cabinet
Model No. 1824	Flammables Cabinet
Model No. 1830	Flammables Cabinet
Model No. 1836	Flammables Cabinet
Model No. 1848	Flammables Cabinet
Model No. 2036	Partially Lined Combination Cabinet
Model No. 2048	Partially Lined Combination Cabinet
Model No. 2236	Fully Lined Combination Cabinet
Model No. 2248	Fully Lined Combination Cabinet



# Under-the-Counter Cabinet

## Installation and Maintenance Instructions

Please unpack these cabinets carefully. If any damage is evident, immediately notify the delivering freight carrier and the supplier of your SciMatCo Under-the-Counter cabinet.

You have just purchased the finest cabinet designed solely for the safe storage of hazardous chemicals. Please read the manual thoroughly before using this product. Only store chemicals in the compartments which are designed to hold these chemicals. Please keep this manual on hand for future reference.

**Caution:** Always use eye, face, hand, and body protection when using this cabinet.

### I. Unpacking the Cabinets

The best method of unpacking the cabinets is as follows.

- A) Remove staples, tape and banding material from the top and bottom of the carton.
- B) Open the top of the box and remove the four (4) packing corners. Gently lift the carton off the product.
- C) Move the cabinet to its appropriate location. Make sure the surface on which you would like to set the cabinet can support the weight of the cabinet and its contents.

You must level the cabinet to insure proper door operation. Four leveling feet located at the bottom of the cabinet are to be used for this purpose.

**Caution:** Do not locate your cabinet near open flame, spark or heat or any source of ignition.

### II. Assembly

Your Under-the-Counter cabinets come fully assembled so no assembly is needed.

**Warning:** Always close and secure cabinet doors after removing your chemicals.

### III. Maintenance

#### Acid Storage Compartment

Your SciMatCo Under-the-Counter acid cabinet is constructed of 1" thick, high density plywood, painted with several coats of a lacquer paint and then painted with three (3) coats of clear epoxy. This cabinet was designed to give you years of trouble free service with minimal maintenance.

**Caution:** Always use eye, face, hand and body protection when cleaning this cabinet.

The inside of the acid cabinet should occasionally be wiped off using a weak solution of sodium phosphate tribasic. Having a good cleaning program and following it will help you get many years of service out of you cabinet. The ratio is 45 grams of sodium phosphate tribasic to one gallon of water.


The cabinets should also be completely emptied three (3) times a year and cleaned thoroughly using this solution. Remember, the greatest stress is placed on the cabinet during the hot summer months. Acid bottles are designed to vent and will vent excessively during periods of high temperatures. Keep this in mind when planning your cleaning schedule. Having a good cleaning program and following it will assure you that you will get many years of service out of your cabinet.

Any scratches in the painted surface of the cabinet should be repainted at once. You can order an epoxy touch up paint kit from us. Please refer to the replacement parts list at the end of this manual for the appropriate catalog number.

The wooden hinges, on your acid cabinet, is designed for long life, however, they may occasionally squeak. If this occurs, a lubricant, such as a petroleum spray can be used on the hinges to relieve this problem.

**Danger:** Do not store perchloric acid in this cabinet or in the company of any other chemicals. Always store perchloric acid alone. Perchloric acid forms explosive crystals, which can cause an explosion.

**Warning:** Absolutely no flammable materials should ever be stored in the acid section of the cabinet. This provision excludes all acids which may also be flammable.



Nitric Acid Compartment as shown in a standard style acid cabinet

## Flammables Compartment

Your Under-the-Counter flammables cabinet is constructed of 1" thick, high density plywood, painted with several coats of a lacquer paint and then painted with three coats of a clear epoxy. This cabinet was designed to give you years of trouble free service with minimal maintenance.

The flammables compartment should be occasionally wiped off with a wet cloth, using warm water and a common household cleaner. Any scratch in the cabinet should also be repaired immediately by touching them up with the epoxy touch up kit which is available from us.

Your flammables compartment has an all metal hinge. From time to time it may squeak. You may use a petroleum base lubricant on this door hinge to relieve this problem.

**Warning:** Absolutely no acid or corrosive materials should ever be stored in the flammables section of the cabinet.

## IV. Nitric Acid Storage (Acid Compartment only)

There is an urgent need to keep the nitric acid separate from the acetic acid in your acid cabinets. The combination of these two substances has resulted in fires and in small explosions. We offer an accessory for our acid cabinets which has been designed especially for the storage of nitric acid. The use of this welded polypropylene container dramatically reduces the possibility of the mixtures of these two dangerous substances. Please use part number SC8071 when ordering and contact your dealer for pricing.

**Warning:** If you do not wish to purchase the optional Nitric Acid Compartment, please do not store the acetic and nitric acid in the same cabinet. A fire or explosion may result.

## V. Replacement Parts List

<b>Description</b>	<b>Part No.</b>
<b>Paint Touch-up Kits</b>	
Grey	2060-1
Blue	8041-7
Yellow	7021-5
<b>Extra Keys</b>	
Acid	8041-4
Flammables/Poison	7021-1
<b>New Locks</b>	
Acid	8041-5
Flammables/Poison	7021-4
<b>Miscellaneous Parts</b>	
TSP Cabinet Cleaner for Acid Cabinet	8041-8

If you have any questions or need any replacement parts, please contact your supplier.

## SciMatCo Product Warranty

Scientific Materials Company warrants its products to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. This warranty excludes damage caused by abuse, misuse, accidents or the failure to provide proper maintenance.

**SciMatCo**  
Division of Flinn Scientific, Inc.  
P.O. Box 305 Batavia, IL 60510  
(888) 879-7620 FAX (630) 879-8893



# Acid/Corrosives Cabinet

## Installation and Maintenance Instructions

Model No. 8041  
Model No. 8051  
Model No. 8063  
Model No. 8081

Partially Lined Model  
Fully Lined Model  
Bench Model  
Econ-A-Cab Model

Please unpack these cabinets carefully. If any damage is evident, immediately notify the delivering freight carrier and the supplier of your SciMatCo cabinet.

You have just purchased the finest cabinet designed solely for the safe storage of corrosive chemicals. Please read the manual thoroughly before using this product. Only store corrosive materials in this cabinet. Please keep this manual on hand for future reference.

**Caution:** Always use eye, face, hand, and body protection when using this cabinet.

### I. Unpacking the Cabinets

The best method of unpacking the cabinets is as follows.

- A) Remove staples, tape and banding material from the top and bottom of the carton.
- B) Open the top of the box and remove the four (4) packing corners. Gently lift the carton off the product. Be careful not to scratch the cabinet with the box.
- C) Move the cabinet to its appropriate location. Make sure the surface on which you would like to set the cabinet can support the weight of the cabinet and its contents.

**Caution:** Do not locate your cabinet near open flame, spark, heat or any source of ignition.

### II. Assembly and Set-Up

Your Acid/Corrosives cabinet comes fully assembled so no assembly is needed. Once you have moved your cabinet into its proper location, you may need to adjust the leveling feet to insure the proper stability.

**Warning:** Always close and secure cabinet doors after removing your chemicals.

### III. Maintenance

Your SciMatCo Acid/Corrosives cabinet is constructed of 1" thick, high density plywood, painted with several coats of a lacquer paint and then painted with three (3) coats of clear epoxy. This cabinet was designed to give you years of trouble free service with minimal maintenance.

**Caution:** Always use eye, face, hand and body protection when cleaning this cabinet.


The inside of the acid cabinet should occasionally be wiped off using a weak solution of sodium phosphate tribasic. Sodium phosphate cleaner can be found at the hardware store. The ratio is 45 grams of sodium phosphate tribasic to one gallon of water.

The cabinets should also be completely emptied three (3) times a year and cleaned thoroughly using this solution. Remember, the greatest stress is placed on the cabinet during the hot summer months. Acid bottles are designed to vent and will vent excessively during periods of high temperatures. Keep this in mind when planning your cleaning schedule. Having a good cleaning program and following it will assure you that you will get many years of service out of your cabinet.

Any scratches in the painted surface of the cabinet should be repainted at once. You can order an epoxy touch up paint kit from us. Please refer to the replacement parts list at the end of this manual for the appropriate catalog number.

The wooden hinges are designed for long life, however, they may occasionally squeak. If this occurs, a lubricant, such as a petroleum spray can be used on the hinges to relieve this problem.

**Danger:** Do not store perchloric acid in this cabinet or in the company of any other chemicals. Always store perchloric acid alone. Perchloric acid forms explosive crystals, which can cause an explosion.



Nitric Acid Compartment may fit on top shelf depending on cabinet.

Nitric Acid Compartment as shown in a standard style acid cabinet

## IV. Nitric Acid Storage

There is an urgent need to keep the nitric acid separate from the acetic acid in your acid cabinet. The combination of these two substances has resulted in fires and small explosions.

We offer an accessory for our acid cabinets (except Model No. 8063) which has been designed especially for the storage of nitric acid. The use of this welded polypropylene container dramatically reduces the possibility of the mixture of these two dangerous substances. For ordering information, please see the replacement parts list.

**Warning:** If you do not wish to purchase the optional Nitric Acid Compartment, please do not store the acetic and nitric acid in the same cabinet. A fire or explosion may result.

## V. Materials to be stored

This Acid/Corrosives cabinet is designed for the storage of corrosive materials, primarily acids and bases.

**Warning:** Absolutely no flammable materials should ever be stored in this cabinet.

By design this cabinet does not have the metal hinges and lock which are required for use as a flammables cabinet. This cabinet has been designed to give years of service as an Acid/Corrosives cabinet by exempting the use of metal in the construction of the cabinet. Never store flammable materials in this cabinet unless they are also corrosive materials.

## VI. Replacement Part List

<b>Description</b>	<b>Part No.</b>
<b>Paint Touch-up Kit</b> Blue	8041-7
<b>Extra Keys</b> Acid/Corrosive	8041-4
<b>New Locks</b> Acid/Corrosive	8041-5
<b>Miscellaneous Parts</b> Dispensing Tray	8041-6
<b>Miscellaneous Parts</b> Nitric Acid Compartment for all cabinets except SC8063	8071

If you have any questions or need any replacement parts, please contact your supplier.

## SciMatCo Product Warranty

Scientific Materials Company warrants its products to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. This warranty excludes damage caused by abuse, misuse, accidents or the failure to provide proper maintenance.

# SciMatCo

Division of Flinn Scientific, Inc.  
P.O. Box 305 Batavia, IL 60510  
(630) 879-7620 (888) 879-7620  
FAX (630) 879-8893  
scimatco@scimatco.com



# Flammables Cabinet

## Installation and Maintenance Instructions

Please unpack these cabinets carefully. If any damage is evident, immediately notify the delivering freight carrier and the supplier of your SciMatCo cabinet.

You have just purchased the finest cabinet designed solely for the safe storage of hazardous chemicals. Please read the manual thoroughly before using this product. Only store chemicals in the compartments which are designed to hold these chemicals. Please keep this manual on hand for future reference.

**Caution:** Always use eye, face, hand, and body protection when using this cabinet.

### I. Unpacking the Cabinets

The best method of unpacking the cabinets is as follows.

- A) Remove staples, tape and banding material from the top and bottom of the carton.
- B) Open the top of the box and remove the four (4) packing corners. Gently lift the carton off the product. Be careful not to scratch the cabinet with the box.
- C) Move the cabinet to its appropriate location. Make sure the surface on which you would like to set the cabinet can support the weight of the cabinet and its contents.

**Caution:** Do not locate your cabinet near open flame, spark or heat or any source of ignition.

### II. Assembly and Set-Up

Your SciMatCo Flammables cabinet comes fully assembled so no assembly is needed. Once you have moved your cabinet into its proper location, you may need to adjust the leveling feet on floor models to insure the proper stability.

**Warning:** Always close and secure cabinet doors after removing your chemicals.

### III. Maintenance

Your SciMatCo Flammables cabinet is constructed of 1" thick, high density plywood, painted with several coats of a lacquer paint and then painted with three (3) coats of clear epoxy. This cabinet was designed to give you years of trouble free service with minimal maintenance.

**Caution:** Always use eye, face, hand and body protection when cleaning this cabinet.

The inside of the cabinet should occasionally be wiped out using a solution of any household cleaner. Having a good cleaning program and following it will assure you that you will get many years of service out of your cabinet.

Any scratches in the painted surface of the cabinet should be repainted at once. You can order an epoxy touch up paint kit from us. Please refer to the replacement parts list at the end of this manual for the appropriate catalog number.



## IV. Materials to be stored

This SciMatCo Flammables cabinet is designed for the storage of flammable liquids.

**Warning:** Absolutely no corrosive materials should ever be stored in this cabinet.

The metal hinges and lock used in this cabinet, while required for use as a flammables cabinet, will not hold up to corrosive materials such as acids. Never store acids or any material which will corrode the metal components in this cabinet.

## V. Venting Your Cabinet

Local Fire Departments may require the venting of Flammables Cabinets. The NFPA, (National Fire Protection Association) which has developed many of the fire codes throughout the United States, sees no benefit in the ventilation of a flammables cabinet for the purpose of fire prevention. If you must conform to a local standard requiring the venting of your cabinet, you must follow these guidelines in order to prevent compromising the fire protection offered by your flammables cabinet.


1) The cabinet must have the proper flame arrestors installed. (Contact SciMatCo for the Flame Arrestors)

2) Vents should be located on opposite sides of the cabinet, one at the top and one at the bottom. The vent at the top should be used as the inlet port and the bottom vent should be used as the exhaust port. (See figure below)

3) Exhaust piping must be rigid steel piping and be connected to the lower vent. If you use a mechanical blower it must be an explosion proof motor blower assembly.

4) Maximum exhaust duct length of 25 feet must be observed.

Failure to observe these guidelines may seriously reduce the fire protection of any flammables cabinet.



## VI. Replacement Parts List

<b>Description</b>	<b>Part No.</b>
<b>Paint Touch-up Kit</b>	
Yellow	7021-5
<b>Flame Arrestors</b>	7021-8
<b>Extra Keys</b>	
Bench, 18 Gallon and Wall mount models	7021-1
45 Gallon Model	7131-1
<b>New Locks</b>	
Bench, 18 Gallon and Wall mount models	7021-4
45 Gallon Model	7131-2

If you have any questions or need any replacement parts, please contact your supplier.

## SciMatCo Product Warranty

Scientific Materials Company warrants its products to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. This warranty excludes damage caused by abuse, misuse, accidents or the failure to provide proper maintenance.

**SciMatCo**

Division of Flinn Scientific, Inc.  
P.O. Box 305 Batavia, IL 60510  
(630) 879-7620 FAX (630) 879-8893



# Poisons Cabinet

## Installation and Maintenance Instructions

Model No. 6080

Bench Model

Please unpack these cabinets carefully. If any damage is evident, immediately notify the delivering freight carrier and the supplier of your SciMatCo cabinet.

You have just purchased the finest cabinet designed solely for the safe storage of hazardous chemicals. Please read the manual thoroughly before using this product. Only store chemicals in the compartments which are designed to hold these chemicals. Please keep this manual on hand for future reference.

**Caution:** Always use eye, face, hand, and body protection when using this cabinet.

### I. Unpacking the Cabinets

The best method of unpacking the cabinets is as follows.

- A) Remove staples, tape and banding material from the top and bottom of the carton.
- B) Open the top of the box and remove the four (4) packing corners. Gently lift the carton off the product. Be careful not to scratch the cabinet with the box.
- C) Move the cabinet to its appropriate location. Make sure the surface on which you would like to set the cabinet can support the weight of the cabinet and its contents. Place cabinet on a surface no higher than 30" from the floor.

### II. Assembly and Set-Up

Your Poisons cabinet comes fully assembled so no assembly is needed.

**Warning:** Always close and secure cabinet doors after removing your chemicals.

### III. Maintenance

Your SciMatCo Poisons cabinet is constructed of 1" thick, high density plywood, painted with several coats of a lacquer paint and then painted with three (3) coats of clear epoxy. This cabinet was designed to give you years of trouble free service with minimal maintenance.

**Caution:** Always use eye, face, hand and body protection when cleaning this cabinet.

The inside of the cabinet should occasionally be wiped out using a weak solution of any household cleaner. Having a good cleaning program and following it will assure you that you will get many years of service out of your cabinet.

Any scratches in the painted surface of the cabinet should be repainted at once. You can order an epoxy touch up paint kit from us. Please refer to the replacement parts list at the end of this manual for the appropriate catalog number.

## IV. Materials to be stored:

This Poisons cabinet is designed for the storage of poisons.

**Warning:** *Absolutely no corrosive materials should ever be stored in this cabinet.*

The metal hinges and lock used in this cabinet, while required for use as a poisons cabinet, will not hold up to strong oxidizers such as calcium hypochlorite, iodine or corrosive materials such as acids.

Never store acids or any material which will corrode the metal components in this cabinet.

## V. Inventory Cards

A package of 50 blank perpetual inventory cards has been included with your cabinet. Each time you elect to store a toxic substance in this cabinet, we urge you to complete one of the preprinted, perpetual inventory cards and file that card in the rack provided within the cabinet.

**Warning:** *The disposal of substances used in connection with this equipment may be governed by various federal, state or local regulations. All users of this equipment are urged to become familiar with any regulations that apply in the user's area concerning the dumping of waste materials in or upon water, land or air and to comply with such regulations.*

## VI. Replacement Parts List

Description	Part No.
<b>Paint Touch-up Kit</b>	
Grey	2060-1
<b>Extra Keys</b>	
Poisons	7021-1
<b>New Locks and Keys</b>	
Poisons	7021-4
<b>Miscellaneous Parts</b>	
Package of inventory cards	6080-2

If you have any questions or need any replacement parts, please contact your supplier:

## SciMatCo Product Warranty

Scientific Materials Company warrants its products to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. This warranty excludes damage caused by abuse, misuse, accidents or the failure to provide proper maintenance.

# SciMatCo

Division of Flinn Scientific, Inc.  
P.O. Box 305 Batavia, IL 60510  
(630) 879-7620 FAX (630) 879-8893