

date:

23.10.2009

Operating Instructions

Electroporator Eppendorf 2510



Hazards for Humans and Environment



Risk of liberation of genetically modified organisms (GMO)
Risk of electrical shock

Protective Measures and Rules of Behaviour

Prepare electrocompetent cells (*i.e.* cells in a salt-free solution): Bacteria, yeast or cell lines

Electrocompetent bacteria or yeasts can be stored at - 20 °C (recommended - 80 °C)

All working steps should be done in an ice-bath:

- Resuspend the vector to be inserted into the cells in 5 μl A. dest. and transfer into a sterile Eppendorf tube
- Pre-cool the electroporation cuvettes: 10 mm gap for bacteria and yeast, 40 mm gap for cell lines
- Transfer sterile for each electroporation 1 ml of medium (LB, YC/Man-ol) into a sterile test tube
- Turn on electroporator, adjust voltage for the pulse: *E. coli* 1,250 V/10 mm, *Sc.c.* 750 V/10 mm
- Thawed electrocompetent micro-organism slowly on ice
- Add 40 μl cells to the 5 μl vector, mix softly by pipetting and transfer all the mixture into the electr. cuvette
- Transfer the electr. cuvette into the electroporator and pulse immediately by pressing twice "PULSE"
- Take off electr. cuvette and add immediately 1 ml of medium (better fast than sterile)
the medium is w/o selection pressure
- Transfer all cells into the Eppendorf tube and incubate at appropriate temp (37 °C or 30 °C) for 1 or 3 h:
transfected cells start expressing the resistance gene introduced with the vector
- Plate 10 μl (vector) or 100 μl (ligation) of the cells onto selective agar plate
- Centrifuge the cells (30 s, 10,000 g) and plate remaining cells onto selective agar plates

Cleaning & Maintenance

Cleaning electroporation cuvette for reuse:

Add 1 ml 70 % EtOH to kill remaining GMO (overnight)

Add 1 ml 1 M HAc to delete residual vector DNA (4 h)

Wash thoroughly (10 x) with A. dest. to eliminate all acid

Add 1 ml 70 % EtOH and store slightly taped

until all EtOH has evaporated

tap pos:

1

1

2



